



nivia®



AN ISO 18001 CERTIFIED COMPANY

A test is worth thousand expert opinions



Instruments for
GENERAL SURGERY
INTESTINAL SURGERY
UROLOGY
GYNECOLOGY



Rajiv Munshi
Managing Director

Aishwariya Meditech Pvt. Ltd.

7/9, Loni Industrial Area, Site-II,
Mohan Nagar, Ghaziabad-201007

- Delhi
- Bhopal
- Bangalore
- Kolkata

+91 120-2983137

+91 9311084747

aishwariya.meditech@gmail.com

www.aishmed.com

• MOT • MGPS • CSSD



Also available, on request
Instruments for
DENTAL SURGERY
ORTHOPAEDIC SURGERY
and **EYE SURGERY**

A Synonym of Quality, Precision & Reliability
Since 1961

nivia[®]

AN ISO 18001 CERTIFIED COMPANY

Our Vision

To become a vibrant country wide acknowledged top league player in hospital furniture manufacturing and allied business, with firm commitment for growth, innovations, quality & customer focus.

Our Mission

To become leading manufacturer of hospital furniture and related businesses by ushering in cost effective, high quality, innovative products and services, enhancing value to esteemed customers, fostering mutually beneficial partnerships during this continuous journey.

Our Values

Building trust with integrity and transparency in all transactions.
Developing on the valuable suggestions from our customers to exceed their expectations.
Continuous R & D to provide innovative, friendly solutions and products,
Surpassing the established benchmarks and defining ever excelling standards for the industry to follow.

The Company

Nivia Industries established in 1961 have been pioneer in central India for manufacturing medical equipments, hospital equipments / furniture and surgical equipments.

The assured quality oriented cost effective products have gained us our customer trust.

Customer's Trust.

We feel proud to be part of an industry that has been serving the health care sector for decades.

Our manufacturing facilities are located at industrial township of Bhopal.

The ever expanding list of our esteemed clients and customers speak volumes of our sincere efforts towards our vision to be in top league in surgical/hospital appliances and allied products.

Our Team

The driving force of Nivia Industries is Mr Sanjeev Sahani, who has hands on experience in manufacturing, and trading of four decades in the field of surgical equipments and hospital furniture.

He has been providing consultancy to leading hospitals.

With his immense experience, clients have been benefitted by his cost effective solutions.

His single minded devotion to work, vast knowledge and insight has been instrumental in galvanizing organization in to a leading one with tremendous growth.

Volume and Diversity.

Nivia Industries is well equipped with latest technology, tools and machinery with highly motivated work force.

Man and machine synchronization is optimized to produce quality products in repetitively.

SURGICAL INSTRUMENTS



**B.P. HANDLE
STANDARD SOLID**
NB-001 No.03



**B.P. HANDLE
STANDARD SOLID**
NB-002 No.04



**B.P. HANDLE STANDARD
SOLID (DOUBLE ENDED)**
NB-003 No.03, 04



TONGUE DEPRESSOR BUCHWALD LACK'S	
STRAIGHT	CURVED
NB-004 5"	NB-005 5"

SURGICAL INSTRUMENTS



**SKIN GRAFTING HANDLE
WATSON MOD. SS**
NB-006



**ARTERY FORCEPS
SPANCER WELLS**

	STRAIGHT	CURVED
NB-007	6"	NB-010 6"
NB-008	7"	NB-011 7"
NB-009	8"	NB-012 8"



**MOSQUITO BABY FORCEPS
HARTMANN**

	STRAIGHT	CURVED
NB-013	4"	NB-014 4"



**MOSQUITO HALSTED
FORCEPS**

	STRAIGHT	CURVED
NB-015	5"	NB-016 5"

SURGICAL INSTRUMENTS



ARTERY FORCEPS (MIXTER) KELLY			
STRAIGHT		CURVED	
NB-017	6"	NB-018	6"

ARTERY FORCEPS (MIXTER) DANDY			
STRAIGHT		CURVED	
N-019	5½"	N-020	5½"



ARTERY FORCEPS (MIXTER) BABY	
N-021	5"

ARTERY FORCEPS MIXTER (CURVED)	
N-022	6"
N-023	8"
N-024	10"

SURGICAL INSTRUMENTS



KOCHER FORCEPS OCHSNER			
STRAIGHT		CURVED	
NB-025	6"	NB-027	6"
NB-026	8"	NB-028	8"



ALLIS TISSUE FORCEPS BOYS	
NB-029	6"



ALLIS TISSUE FORCEPS	
NB-030	6"
NB-031	8"



ALLIS TISSUE FORCEPS DEBAKEY	
NB-032	6"
NB-033	8"



SURGICAL INSTRUMENTS



BABCOCK FORCEPS BABY

NB-034	5"
--------	----



BABCOCK FORCEPS

NB-035	6"
NB-036	8"



**GLAND HOLDING FORCEPS
LANE**

NB-037	6"
NB-037	7"
NB-037	8"



**NEEDLE HOLDER
CRILE WOOD**

NB-038	4"
NB-039	6"
NB-040	8"

SURGICAL INSTRUMENTS



NEEDLE HOLDER MAYO HEGAR	
NB-041	6"
NB-042	7"
NB-043	8"



NEEDLE HOLDER MASSON	
NB-044	9"



NEEDLE HOLDER BOZEMANN	
NB-045	5"
NB-046	6"



NEEDLE HOLDER WITH SCISSORS OLSEN-HEGAR	
NB-047	5½"

SURGICAL INSTRUMENTS



LIGATURE NEEDLE SYME
NB-048 6½"



LIGATURE ADJUSTABLE CASPAR
NB-049 9"



PROBE MALLEABLE			
	SS	COPPER	
NB-050	7"	NB-051	7"
NB-050	8"	NB-051	8"
NB-050	10"	NB-051	10"



PROBE WITH EYE	
NB-052	6"
NB-052	7"
NB-052	8"

SURGICAL INSTRUMENTS



**PROBE (FISTULA)
LOCKHART-MUMMERY**

NB-053	6 1/2"
--------	--------



**SCISSORS FINE
SHARP & BLUNT**

	STRAIGHT	CURVED
NB-054	6"	NB-056 6"
NB-055	8"	NB-057 8"



SCISSORS METZENBAUM

	STRAIGHT	CURVED
NB-058	4"	NB-062 4"
NB-059	6"	NB-063 6"
NB-060	8"	NB-064 8"
NB-061	9"	NB-065 9"



**SCISSORS METZENBAUM
FINE**

	STRAIGHT	CURVED
NB-066	6"	NB-068 6"
NB-067	8"	NB-069 8"

SURGICAL INSTRUMENTS



SCISSORS MAYO

STRAIGHT		CURVED	
NB-070	6"	NB-072	6"
NB-071	8"	NB-073	8"



SCISSORS KILNER

STRAIGHT		CURVED	
NB-074	4 1/4"	NB-076	4 1/4"
NB-075	6"	NB-077	6"



SCISSORS JOSEPH

STRAIGHT		CURVED	
NB-078	5 1/2"	NB-079	5 1/2"



SCISSORS BONN

STRAIGHT		CURVED	
NB-080	3 1/2"	NB-082	3 1/2"
NB-081	4 1/2"	NB-083	4 1/2"

SURGICAL INSTRUMENTS



SCISSORS LIGATURE SPENCER
 NB-084 5½"



SCISSORS STITCH CUTTING
 NB-085 6"



SCISSORS REYNOLDS CURVED
 NB-086 5"
 NB-087 6"



SCISSORS UMBILICAL MOD. USA
 NB-088 4¼"



SURGICAL INSTRUMENTS



**SCISSORS
BANDAGE CUTTING**
NB-089 7"



**SCISSORS EPISIOTOMY
BRAUN-STADLER**
NB-090 5 1/4"



**DISSECTING FORCEPS
STANDARD SERRATED**
NB-091 6"
NB-092 8"



**DISSECTING FORCEPS
SERRATED POTTS-SMITH**
NB-093 7"

SURGICAL INSTRUMENTS



**DISSECTING FORCEPS
STANDARD 1:2 TEETH**

NB-094	6"
NB-095	8"



**ORGAN HOLDING FORCEPS
DUVAL**

NB-096	7"
	OPTIONAL 8"



**DISSECTING FORCEPS
ADSON SERRATED**

NB-097	4"
NB-098	5"



**DISSECTING FORCEPS
ADSON 1:2 TEETH**

NB-099	4"
NB-100	5"

SURGICAL INSTRUMENTS



**DISSECTING FORCEPS
RUSS. MODELL**

NB-101	6"
NB-102	8"



**INTESTINAL FORCEPS
DENIS-BROWNE**

NB-103	7½"
--------	-----



**TISSUE FORCEPS
GRAEFE MULTI TOOTH**

NB-104	4½"
--------	-----



**RETRACTOR SKIN
SINGLE HOOK VICKERS**

NB-105	7"
--------	----

SURGICAL INSTRUMENTS



**RETRACTOR SKIN
TWO HOOK VICKERS**
NB-106 7"



**RETRACTOR SKIN
THREE HOOK VICKERS**
NB-107 7"



**RETRACTOR SKIN
GILLIES**
NB-108 8"



**RETRACTOR SKIN
SENN-MUELLER**

	SHARP	
NB-109		6 1/4"
	BLUNT	
NB-110		6 1/4"

SURGICAL INSTRUMENTS



RETRACTOR VOLKMANN

SHARP	8 1/4"
NB-111	1 PRONG
NB-112	2 PRONGS
NB-113	3 PRONGS
NB-114	4 PRONGS
BLUNT	8 1/4"
NB-115	1 PRONG
NB-116	2 PRONGS
NB-117	3 PRONGS
NB-118	4 PRONGS



RETRACTOR DOYEN WITH RING HANDLE (8 1/2")

NB-119	2"
NB-120	2 1/2"
NB-121	3"
NB-122	3 1/2"
NB-123	4"

RETRACTOR LANGENBECK 8 1/2"

NB-124	6 X 25 mm
NB-125	8 X 35 mm
NB-126	11 X 35 mm
NB-127	15 X 35 mm
NB-128	11 X 41 mm
NB-129	11 X 55 mm
NB-130	14 X 70 mm
NB-131	12 X 80 mm
NB-132	16 X 80 mm



SURGICAL INSTRUMENTS



**RETRACTOR MORISS
WITH RING HANDLE (8 1/2")**

NB-133	2"
NB-134	2 1/2"
NB-135	3"
NB-136	3 1/2"
NB-137	4"



**RETRACTOR MINI LAP
DEAVER**

NB-138	No. 1
NB-139	No. 2
NB-140	No. 3



RETRACTOR DEAVER

NB-141	1"
NB-142	1 1/2"
NB-143	2"
NB-144	2 1/2"
NB-145	3"



RETRACTOR CZERNY

NB-146	6"
NB-147	7"
NB-148	8"

SURGICAL INSTRUMENTS



RETRACTOR ROUX

NB-149	(24 x 20/20 x 28mm)	15cm/6"
NB-150	(29 x 26/36 x 36mm)	6cm/6 1/4"
NB-151	(33 x 29/38 x 43mm)	7cm/6 3/4"



RETRACTOR LANGENBECK

N-152	8 1/2" - 60X 20 mm
-------	--------------------



RETRACTOR IMONS

NB-153	115 X 22 mm
NB-154	115 X 27 mm



RETRACTOR KELLY

NB-155	160 X 38 mm
NB-156	160 X 57 mm

SURGICAL INSTRUMENTS



**SELF RETAINING
RETRACTOR PLESTER 2X2**

NB-157	5" SHARP
NB-158	5" BLUNT



**SELF RETAINING
RETRACTOR WEITLANER 2X3**

NB-159	4 1/4" SHARP
NB-160	4 1/4" BLUNT



**SELF RETAINING
RETRACTOR MOLLISON 4X4**

NB-161	6" SHARP
--------	----------



**SELF RETAINING
RETRACTOR BALFOUR BABY**

NB-162	10" SHAFT
--------	-----------



SURGICAL INSTRUMENTS



**SELF RETAINING RETRACTOR
BALFOUR STANDARD**

NB-163	8" SHAFT
NB-164	10" SHAFT



**SELF RETAINING RETRACTOR
BALFOUR MOD.**

N-165	10" SHAFT
-------	-----------



**SELF RETAINING RETRACTOR
BALFOUR STANDARD**

NB-166	8" SHAFT
NB-167	10" SHAFT



**SELF RETAINING RETRACTOR
BALFOUR MOD.
(DEEP LATERAL BLADE)**

NB - 168	10" SHAFT
----------	-----------

SURGICAL INSTRUMENTS



**SELF RETAINING RETRACTOR
POZZI'S**
NB-169




**INTESTINAL CLAMP
BABY KOCHER**

STRAIGHT	CURVED
NB-170 5"	NB-171 5"



**INTESTINAL CLAMP
DOYEN**

STRAIGHT	CURVED
NB-172 10"	NB-173 10"



PROCTOSCOPE KELLY

NB-174	SMALL
NB-175	MEDIUM
NB-176	LARGE
NB-177	EXTRA LARGE

SURGICAL INSTRUMENTS



RECTAL SPECULA
KILLIAN (FOR CHILDREN)

NB-178

5½"



HAEMORRHOIDAL LIGATOR
WITH LOADING ATTACH. &
SUCTION PORT

NB -179



SUCTION CANULA
YASARGIL (CP)

NB-180

1,2,3,4



HYSTERECTOMY CLAMP
WERTHEIM-CULLEN

NB-181

8½"

SURGICAL INSTRUMENTS



**SCISSORS POTTS-DE MARTELL
GALL DUCT**

NB-182 9"



RECTAL DILATOR ST. MARKS

NB-183	No. 1
NB-184	No. 2 or 3



**PROBE GALL STONE
MAYO MALLEABLE**

NB-185 13"



**DILATOR BAKES MOD.
MALLEABLE
(COMMON BILE DUCT DILATOR)**

NB-186 12"

SURGICAL INSTRUMENTS



6/10FG, 8/12FG, 10/14FG, 12/16FG,
14/18FG, 16/20FG, 18/22FG, 20/24FG,
22/26FG, 24/28FG, 26/30FG, 28/32 FG

**DILATOR URETHRAL
CLUTTON SET OF 12**
NB-187 10½"



1/4FG, 2/5FG, 3/6FG, 4/7FG, 5/8FG,
6/9FG, 7/10FG, 8/11FG, 9/12FG,
10/13FG, 11/14FG, 12/15FG

**DILATOR URETHRAL
LISTER WITH OLIVE ENDS
SET OF 12**
NB-188 11"



**KIDNEY ELEVATING FORCEPS
DARTIGUES**
NB-189 10"



**KIDNEY ELEVATING FORCEPS
SOMER**
NB-190 10"

SURGICAL INSTRUMENTS



KIDNEY STONE FORCEPS RANDALL

NB-191	23cm/9" (No. 1)
NB-192	23cm/9" (No. 2)
NB-193	23cm/9" (No. 3)
NB-194	23cm/9" (No. 4)



KIDNEY STONE FORCEPS DESJARDINS

NB-195	No. 1
NB-196	No. 2



URETER FORCEPS DEBAKEY

NB-197	9 1/2"
--------	--------

SURGICAL INSTRUMENTS



CATHETER

MALE		FEMALE	
NB-198	No. 1	NB-201	No. 1
NB-199	No. 2	NB-202	No. 2
NB-200	No. 3	NB-203	No. 3



INTRODUCER SLOTTED (FOR FOLEY)

NB-204



TROCAR NELSON MOD. SLOTTED

NB-205 No. 1, 2, 3, 4



BULLDOG CLAMP DE-BAKEY

STRAIGHT		CURVED	
NB-206	3"	NB-208	3"
NB-207	3½"	NB-209	3½"

SURGICAL INSTRUMENTS



RASPATORIES DOYEN

CHILDREN 7"		ADULT 7"	
NB-210	RIGHT	NB-212	RIGHT
NB-211	LEFT	NB-213	LEFT



RIB SPREADERS FINOCCHIETTO

NB-214	INFANT
NB-215	CHILD
NB-216	(MEDIUM) ADULT



SUCTION CANULA POOL STRAIGHT ALL SS

NB-217	8mm
NB-218	10mm



SUCTION CANULA YANKAUER CP

NB-219	8mm
--------	-----



SURGICAL INSTRUMENTS



COLLIN FORCEPS
 NB-220 6 3/4"



MAGILL INTRODUCING FORCEPS

NB-221	INFANT 6"
NB-222	CHILD 8"
NB-223	ADULT 10"



**TOWEL CLIP CROSS ACTION
 ANGLED ENGL. MOD.**
 NB-224 3 3/8"



**TOWEL FORCEPS
 BACKHAUS**

NB-225	4"
NB-226	5"
NB-227	6"
NB-228	7"

SURGICAL INSTRUMENTS



CHEATAL FORCEPS
NB-229 10½"



CHEATAL LH PATTERN FORCEPS
NB-230 11"



**VAGINAL SPECULA
CUSCO STANDARD**

NB-231	SMALL
NB-232	MEDIUM
NB-233	LARGE



**VAGINAL SPECULA
SIMS STANDARD**

NB-234	SMALL
NB-235	MEDIUM
NB-236	LARGE
NB-237	EXTRA LARGE

SURGICAL INSTRUMENTS



**UTERINE CURETTE BLAKE
SHARP & BLUNT**
NB - 238



UTERINE SCOOP GOURDET

NB-239	6mm
NB-240	8mm
NB-241	10mm
NB-242	12mm
NB-243	14mm



BLADDER SOUND COLLIN
NB-244



UTERINE SOUND SIMS
NB-245

SURGICAL INSTRUMENTS



UTERINE DEPRESSOR SIMS
 NB-246 10 3/4"



FLUSHING CURETTE FOR SUCTION
 NB-247 12 1/4"



HYSTERECTOMY CLAMP HEANEY MOD. SINGLE, DOUBLE HOOK

STRAIGHT		CURVED	
NB-248	8"	NB-249	8"



HYSTERECTOMY CLAMP MAINGOT MOD.

STRAIGHT		CURVED	
NB-250	8"	NB-253	8"
NB-251	10"	NB-254	10"
NB-252	12"	NB-255	12"

SURGICAL INSTRUMENTS



**HYSTERECTOMY CLAMP
NIVIA SPECIAL CURVED**
NB-256 10"



**HYSTERECTOMY CLAMP
HEANEY ATRAUMA**
NB-257 9"ST.
NB-258 8"CUR.
NB-259 8"EXTRA CUR.
NB-260 8"ANGLED



**DILATOR UTERINE
HEGAR BLUNT (SET OF 10)**
NB-261 8"
2.5/3mm 7.5/8mm
3.5/4mm 8.5/9mm
4.5/5mm 9.5/10mm
5.5/6mm 10.5/11mm
6.5/7mm 11.5/12mm



**DILATOR VAGINAL
HAWKIN AMBLER SINGLE ENDED**
NB-262-18 18cm/7"
DIA
6 mm 11 mm
7 mm 12 mm
8 mm 13 mm
9 mm 14 mm
10 mm 15 mm

SURGICAL INSTRUMENTS



UTERINE GRASPING FORCEPS
N-263



INTRA UTERINE LEECH- WILKINSON 10 1/4"

NB-264	SMALL
NB-265	MEDIUM
NB-266	LARGE



UTERINE INTRA FRACER 10 1/4"

NB-267	SMALL
NB-268	MEDIUM
NB-269	LARGE



MTP CANNULA SET OF 6
NB-270 10 1/4"

SURGICAL INSTRUMENTS



MTP CANULA SS & CP
10 1/2"
NB-271 4,6,8,10,12



CURRIE CANULA
NB-272 13"



RUBIN CANULA
NB-273 10"



**UTERINE ELEVATOR
MANIPULATOR MOD.
(DOUBLE ACTION)**
NB-274 11"

SURGICAL INSTRUMENTS



**UTERINE ELEVATOR
TEALE-KAHN**
NB-275 11"



**UTERINE TENACULUM
FORCEP BRAUN (SLENDER TYPE)**
NB-276 10"



VULSELLUM FORCEPS TEALE
NB-277 10"



**UTERINE TENACULUM
FORCEP CZERNY**
NB-278 8"

SURGICAL INSTRUMENTS



**UTERINE TENACULUM
FORCEPS SCHRODER**
NB-279 10 3/4"



**UTERINE HOOK
SAUNDERS**
NB-280 13"



**T TYPE HANDLE DOYEN
MYOMA SCREW**
NB-281 7"



**RING TYPE HANDLE DOYEN
MYOMA SCREW**
NB-282 7"

SURGICAL INSTRUMENTS



**CERVICAL BIOPSY FORCEPS
VAN - DOREN**
NB-283 10"



**OVUM FORCEPS
HEYWOOD SMITH 10"**
NB-284 No. 0,00,1,2,3,4



**SPONGE HOLDING
FORCEPS STRAIGHT
FOERSTER**
NB-285 10"



SINUS FORCEPS LISTER
NB-286 6"
NB-287 8"

SURGICAL INSTRUMENTS



**HYSTERECTOMY CLAMP
GWILLIAMS STRAIGHT, CURVED**
NB-288 8"



DIALATORS VASSEL
N-289 7"



GREEN ARMYTAGE FORCEPS
NB-290 8 3/4"



**OBSTETRICS-OBSTETRICAL
FORCEPS SIMPSON BRAUN**
NB-291 14"

SURGICAL INSTRUMENTS



**OBSTETRICS-OBSTETRICAL
FORCEPS ANDERSON**

NB-292 14"



**OBSTETRICS-OBSTETRICAL
FORCEPS KEDARNATH DAS**

N-293 14"



**SCISSORS MAYO
TUNGSTEN CARBIDE**

STRAIGHT		CURVED	
NB-294	6"	NB-296	6"
NB-295	8"	NB-297	8"



**NEEDLE HOLDER MAYO
HEGAR TUNGSTEN CARBIDE**

NB-298	6"
NB-299	7"
NB-300	8"

SURGICAL INSTRUMENTS



**SCISSORS METZENBAUM
TUNGSTEN CARBIDE**

STRAIGHT		CURVED	
NB-301	4"	NB-304	4"
NB-302	6"	NB-305	6"
NB-303	8"	NB-306	8"

DNC SET
NB-307



**F.A. VALVE WITH
FLOW METER**
NB - 308



**HUMIDIFIER BOTTLE
POLYCARBONATE**
NB - 309

SURGICAL INSTRUMENTS



INFANT AMBU BAG SILICONE

AMBU BAG	
NB-310	INFANT
NB-311	CHILD
NB-312	ADULT



CHILD AMBU BAG SILICONE



ADULT AMBU BAG SILICONE



MOX REGULATOR OXYGEN SINGLE GAUGE
NB - 313



MOX REGULATOR OXYGEN DOUBLE GAUGE
NB - 314



MOX REGULATOR NITROUS SINGLE GAUGE
NB - 315

SURGICAL INSTRUMENTS



LARINGOSCOPE-MACINTOSH
ADULT WITH FOUR BLADES

NB - 316

LARINGOSCOPE-MACINTOSH
ADULT WITH THREE BLADES

NB - 317



LARINGOSCOPE-MILLER
PEAD. WITH TWO BLADES

NB - 318

CLEANING & MAINTENANCE OF SURGICAL INSTRUMENTS



nivia manufactures all the surgical instruments under the supervision of specially trained technical staffs and all the instruments go through intensive checking before final packaging. Still we emphasize mostly on cleaning and maintenance to increase the life and working capability of each single instrument.

nivia suggests some interconnected care and maintenance guidelines to all our buyers.

MANUAL CLEANING (MANUAL CLEANING UNDERGOES SOME CAREFUL STAGES)



- ◆ Rinse all the instruments immediately after operation under flowing water to remove all blood, fluids, tissues and other particles.
- ◆ Submerge instruments into a suitable disinfectant with active cleaning properties so that all surfaces, inner cavities, lumens and openings come into contact with the solution.
- ◆ The solution should preferably be a combination of distilled water and neutral PH (7) detergent.
- ◆ Use stiff plastic brushes to clean the remaining and sticky blood stains particularly from hinges, joints, locks, boxes and other malleable parts and working points Wire brushes may be used in a delicate way for grips, handles, within serrations and jaws (not causing serious damages in working points).
- ◆ Examine carefully whether instruments are visibly clean and free from bloodstains and tissues.
- ◆ Test and make sure the instruments are functioning properly.
- ◆ After scrubbing rinse instruments thoroughly under running water. Final rinsing to be done with distilled or deionized water.
- ◆ Dry instrument with an absorbent, soft and lint-free cloth.
- ◆ Dry the instruments in air and store them in a clean and dry environment in instrument box

AUTOCLAVING (AUTOCLAVING IS ESSENTIAL PART BEFORE AND AFTER ALL SURGERIES)



Lubricate all instruments properly.

Steam sterilize using a prevacuum cycle for 4 minutes at a minimum temperature of 132°C(270°F). *Note Where there is a concern about contamination, the World Health Organization amends processing through autoclave cycle for 18 minutes at 134°C(273°F)*

Do not overload the autoclave chamber

Drying time will vary according to load size and should be increased for larger loads

For better life and working capability of each single instrument.

- ◆ Always use the proper sterilization /Cleaning technique to render the instrument in required condition for use.
- ◆ Use ultra sound cleaning machine and drying machine.



FRAGILE



AISHWARIYA MEDITECH PVT LTD

7/9, Loni Industrial Area, Site-II
Mohan Nagar, Ghaziabad U.P. 201 007

www.aishmed.com



- Delhi
- Bhopal
- Bangalore
- Kolkata

Rajiv Munshi
Managing Director

Aishwariya Meditech Pvt. Ltd.

7/9, Loni Industrial Area, Site-II,
Mohan Nagar, Ghaziabad-201007

+91 120-2983137

+91 9311084747

aishwariya.meditech@gmail.com

www.aishmed.com

• MOT • MGPS • CSSD

